

# HOME HEALTHCARE

<del>\* \* \* \* \* \*</del>

LIFE ON with OMRON
Healthcare Monitoring Device

# WARRANTY REGISTRATION

# OMRON

# GET YEA EXTENDED WARRANTY & EARN REWARDS









Scan

Save

Share



Scan Now To Avail The Offer!





No need to secure bills



Additional warranty for 1 year



Completely safe and secure online registration



Hassle-free paperless process



No requirement of warranty card for service















# INSTALL OMRON CONNECT APP TO TRACK YOUR HEALTH

## OMRON

# GIVE YOURSELF THE GIFT OF GOOD HEALTH

GET REWARDED TODAY!

Make your health journey a more fulfilling one by "DOWNLOADING OMRON CONNECT APP TODAY"





#### **HEALTH GIFT**

REWARD YOURSELF FOR EVERY MEASUREMENT YOU TAKE WITH YOUR BLOOD PRESSURE/BODY COMPOSITION MONITOR.



#### BP DIARY

A USEFUL AND
CONVENIENT FUNCTION
FOR DAILY BLOOD
PRESSURE
MANAGENENT.



#### **OMRON CLOUD**

ACCESS YOUR DATA
READILY ON ANY
MOBILE DEVICES BY
LOGGING IN TO YOUR
CLOUD ACCOUNT.

#### TO GET STARTED WITH HEALTHGIFT





















#### **Omron Connect Benefits**

From eliminating manual BP recordings to providing you with insightful data, OMRON connect offers many benefits that will transform the way you manage your heart health.



#### Connected Effortlessly Using Blutooth

Using Bluetooth, transfer measurements to your smartphone seamlessly after a simple, guided set-up process.



#### **Access Anytime, Anywhere**

Get your latest activity and personal health history at your fingertips through the intuitive dashboard. Share this information with your doctor to receive the most appropriate care.



#### **Greater Insight for Better Heart Health Progress**

Heart health data is displayed in clear and insightful graphs to help you better understand your health trends for more effective progress.



#### Record, Manage and Track Health Data Easily

Record your lifestyle habits, gain personal health insights, get medication reminders and create progress reports instantly that can be shared via email all through the BP diary feature.



#### Earn Rewards

HealthGift program rewards you with gifts including Grab or grocery vouchers when you measure your blood pressure regularly.

# OMRON

**OMRON** COMPLETE **HOME-BASED** AUTOMATIC **UPPER ARM** 



# **REDUCE THE CHANCES OF STROKE**





#### **SPECIFICATIONS**

Model HEM-7530T-AP

BP reading display LCD digital display

Transmission method Bluetooth®low energy technology

Wireless communication Frequency range: 2.4 GHz (2400 - 2483.5 MHz) Modulation: GFSK Effective radiated power: <20 dBm

Operation mode Continuous operation

IP classification IP 21

Rating DC6 V 4 W

4"AA" alkaline batteries 1.5 V Power source

**Battery life** Approximately 300 measurements (using new

alkaline batteries)

Cuff circumference 7" to 17" (170 to 420 mm) (included arm cuff: 9" to 17" (220 to 420 mm) applicable to the monitor

Stores up to 90 blood pressure readings

Internal Memory

**Applied part** Type BF (arm cuff), Type CF (electrodes)

**Protection against** Internally powered ME equipment

electric shock

Monitor, Wide-Range D-Ring Cuff, storage case, 4 "AA" Contents alkaline batteries, instruction manual, quick start guide









Gives a clear message if the risk of AF has been detected in patients or not.

OMRON Complete - OMRON'S 1<sup>st</sup> Home Based Heart Health Monitoring

Measure → Record → Share



**FEATURES** 

**BENEFITS** 



Home blood pressure monitor and single lead ECG that gives easy, instant and accurate results



Screens for a wide range of heart conditions (ARRHYTHMIA, NORMAL SINUS RHYTHM, etc.)

Comprehensive screening



Intuitive reading

Intuitive reading: Easy to read results with a clear indicator for normal/abnormal reading.



To access cardiac data, understanding patient health better

A resource for doctor

# BLOOD PRESSURE MONITOR

(UPPER-ARM TYPE)



#### **Product Specifications**

Measurement Range :

Pressure 0 to 299 mmHg Pulse 40 to 180 beats/min

Memory:

Last measurement

Accuracy:

Pressure <u>+</u> 3mmHg Pulse <u>+</u> 5% of display reading

#### HEM - 7121J



#### **Product Specifications**

Measurement Range:

Pressure 0 to 299 mmHg Pulse 40 to 180 beats/min

Memory:

Last measuremen

Accuracy:

Pressure ±3mmHg
Pulse +5% of display reading





#### Value Added Features

Intellisense Technology

Memory of Last Reading

Hypertension Indicator





#### **Value Added Features**

Intellisense Technology Memory of Last Reading Cuff Wrapping Guide







#### Value Added Features

**Body Movement Detection** 

Intellisense Technology

**Hypertension Indicator** 

Irregular Heartbeat detection







CLINICALLY VALIDATED



#### Value Added Features

Bluetooth Connectivity

**Cuff Wrapping Guide** 

Pulse

Intelli™ Sense

**Hypertension Indicator** 

14 Memory





CLINICALLY VALIDATED



#### Value Added Features

**Bluetooth Connectivity** 

**Body Movement Detection** 

**Cuff Wrapping Guide** 

Pulse

Intelli™ Sense

**Hypertension Indicator** 

14 Memory



CLINICALLY VALIDATED









# HEM - 7120-AP INR 2,490/-

#### **Product Specifications**

Measurement Range:

Pressure 0 to 299 mmHg

Memory:

Pressure ±3mmHg Pulse ±5% of display reading

#### HEM-7140T1-AP



#### **Product Specifications**

Pulse 40 to 180 beats/min

Pressure  $\pm 3$ mmHg Pulse  $\pm 5\%$  of display reading

#### HEM-7141T1-AP



#### **Product Specifications**

Pulse 40 to 180 beats/min

Pressure <u>+</u>3mmHg Pulse <u>+</u>5% of display reading

#### **HEM - 7142T1-AP**



#### **Product Specifications**

#### Measurement Range :

Pressure 0 to 299 mmHg Pulse 40 to 180 beats/min

#### Memory:

14 measuremen

#### Accuracy

Pressure  $\pm 3$ mmHg Pulse  $\pm 5\%$  of display reading

#### HEM - 7143T1-AP



#### **Product Specifications**

#### Measurement Range:

Pressure 0 to 299 mmHg Pulse 40 to 180 beats/min

#### Memory:

30 measurement

#### Accuracy:

Pressure  $\pm 3$ mmHg Pulse  $\pm 5\%$  of display reading

#### HEM - 7143T1-AIN



#### **Product Specifications**

#### Measurement Range:

Pressure 0 to 299 mmHg Pulse 40 to 180 beats/min

#### Memory:

30 measurement

#### Accura

Pressure ±3mmHg Pulse ±5% of display reading









Body Movement Indicator

Hypertension Indicator

Irregular Heartbeat detection

Intellisense Technology 14 Memory

**Cuff Wrapping Guide** 

Bluetooth Conectivity

Pulse



CLINICALLY VALIDATED



TECHNOLOGY
BY OMRON HEALTHCARE
IN JAPAN

**Cuff Wrapping Guide** 











Hypertension Indicator

Irregular Heartbeat detection

Intellisense Technology

30 Memory

Bluetooth Conectivity

Pulse



CLINICALLY VALIDATED



TECHNOLOGY BY OMRON HEALTHCARE IN JAPAN











Cuff Wrapping Guide

Body Movement detection

Hypertension Indicator

Irregular Heartbeat detection

Intellisense Technology 30 Memory

30 Memory

Bluetooth Conectivity

Pulse



CLINICALLY VALIDATED



#### Value Added Features

**Cuff Wrapping Guide** 

**Body Movement Indicator** 

**Hypertension Indicator** 

Irregular Heartbeat detection

Average of last 3 readings

Intellisense Technology

360° Accuracy IntelliWrap™ Cuff



CLINICALLY VALIDATED



### Intelli Wrap"





### Value Added Features

**Cuff Wrapping Guide** 

**Body Movement Indicator** 

**Hypertension Indicator** 

Irregular Heartbeat detection

Average of last 3 readings

60 Memory

Intellisense Technology

360° Accuracy IntelliWrap™ Cuff



CLINICALLY VALIDATED





# Onterior









### **Value Added Features**

**Body Movement Indicator** 

**Blood Pressure** 

**Detects Irregular Heartbeat** 

**Date & Time Display** 

Intellisense Technology



#### MOST RECOMMENDED BRAND DOCTORS

CLINICALLY VALIDATED



#### **HEM - 7156-AP**



#### **Product Specifications**

Measurement Range:

Pressure 0 to 299 mmHg Pulse 40 to 180 beats/min

Memory:

60 measurement

Accuracy:

Pressure <u>+</u> 3mmHg Pulse <u>+</u> 5% of display reading

#### **HEM - 7156T-AP**



#### **Product Specifications**

Measurement Range:

Pulse 40 to 180 beats/min

Memory:

60 measurement

Accuracy:

Pressure  $\pm 3$ mmHg Pulse  $\pm 5\%$  of display reading

JPN - 600



INR 5,590/-

#### **Product Specifications**

Measurement Range:

Pressure 0 to 299 mmHg Pulse 40 to 180 beats/min

Memory:

50 measurement

Accuracy :

Pressure ±3mmHg Pulse ±5% of display reading



#### **Product Specifications**

Pulse 40 to 180 beats/min

#### Memory:

Pressure ±3mmHg Pulse ±5% of display reading

#### **HEM - 7600T**



#### **Product Specifications**

Pressure <u>+</u> 3mmHg Pulse <u>+</u> 5% of display reading















#### Value Added Features

Pulse Indicator

**Body Movement Indicator Hypertension Indicator** 

Irregular Heartbeat detection Average of last 3 readings

Intellisense Cuff **AFIB Detection** 

2 Users X 100 Measurements in Memory



CLINICALLY VALIDATED



TECHNOLOGY BY OMRON HEALTHCARE IN JAPAN

# Value Added Features **Body Movement Detection**







**Detects Irregular Heartbeat Cuff Wrapping Guide** Intellisense Technology

360° Accuracy IntelliWarap™ Cuff





**Hypertension Indicator** 




VALIDATED



# Tips For Accurate Measurement

- Measure at the same time of day every day (Within 1 hour after getting up and before bed-time are recommended)
- No bathing, drinking alcohol or caffeine, smoking, exercising or eating 30 minutes before taking a measurement
- \* Use appropriate arm cuff of small, medium, large & wide size

# **Correct Posture For Upper ARM Blood Pressure Measurement**




# **Correct Posture For Wrist Blood Pressure Measurement**

- Place the wrist cuff over wrist.
- · Your palm should face upward.
- Wrap the wrist cuff around your wrist.
- Sit on chair with your feet flat on the floor.
- Sit upright with your back straight.
- The cuff should be at the same level as your heart

Relax your wrist and hand. Do not bend your wrist back, clench your fist, or bend your wrist forward.

Note: make sure that the wrist cuff does not cover the protruding part of the wrist bone (ulna) on the outside of the wrist.



**HEM - 6161** 



#### Product Specifications

#### Measurement Range:

Pressure 0 to 299 mmHa Pulse 40 to 180 beats/min

#### Memory:

Pressure ±3mmHg Pulse ±5% of display reading

#### **HEM - 6181-AP**



#### Product Specifications

Pulse 40 to 180 beats/min

Pressure <u>+</u> 3mmHg Pulse <u>+</u> 5% of display reading

### HEM - 6232T-AP



#### **Product Specifications**

Pulse 40 to 180 beats/min

#### Memory:

Pressure <u>+</u> 3mmHg Pulse <u>+</u> 5% of display reading

# **BLOOD PRESSURE MONITOR** (WRIST TYPE)

#### Value Added Features

Hypertension Indicator

**Detects Irregular Heartbeat** 

Intellisense Technology









TECHNOLOGY BY OMRON HEALTHCARE IN JAPAN

#### Value Added Features

**Body Movement Indicator** 

**Hypertension Indicator** 

**Detects Irregular Heartbeat** 

**Date & Time Display** Morning Average BP

Intellisense Technology







CLINICALLY VALIDATED



TECHNOLOGY BY OMRON HEALTHCARE IN JAPAN

## Value Added Features

**Body Movement Indicator** 

**Hypertension Indicator** 

Detects Irregular Heartbeat

Morning Average BP

Average of last 3 readings

Intellisense Technology

**Dual User** 

Bluetooth Conectivity



CLINICALLY VALIDATED







# **NEBULIZERS**

#### Value Added Features

**High Nebulization rate** 

Effective nebulization

Easy to clean

Optimal particle Size - 4.5 µm





#### Value Added Features

**Effective Nebulization** 

Easy to Use and Clean

**High Nebulization rate** 

Optimal particle Size - 3.9 μm





### Product Specifications

Medication Cup Capacity:

2ml minium and 10 ml maximum

Particle Size:

Nebulizer: 4.56 µm

#### Contonto.

Compressor, Nebulizer Kit, Air tube (PVC, 100cm), Mouthpiece, Adult Mask, Child mask, Air Filer X 2, Instruction Manual



#### **Product Specifications**

Medication Cup Capacity:

Maximum of 12r

Particle Size:

3,90 µm ap

#### Contents:

Compressor, Nebulizer Kit, Air Tube, (PVC, 100cm), Mouthpiece, Nosepiece, Adult Mask (PVC), Child Mask (PVC), 2 pcs Spare Air Filters, Instruction Manual, Quick Start Guide.





#### **Product Specifications**

**Medication Cup Capacity:** 

7 ml (cc) max

Particle Size:

3 µm

Contents:

Compressor, Nebulizer Kit, Air Tube, 5 pcs Replacement Filters, Mouthpiece, Adult Mask, Child Mask, Carry Bag, Instruction Manual



#### **Product Specifications**

Medication Cup Capacity:

1ml minimum and 12 ml maximum

Dantiala Cina

Kit 1:- 4 µm Kit 2:- 9 µm

Contents:

Compressor, Nebulizer Kit, Air tube Mouthpiece, Child mask, Nasal Aspirator, Instruction Manual

#### **NE - U100-AP**



INR 9,790/-

#### **Product Specifications**

**Medication Cup Capacity:** 

Maximum of 10ml

Particle Size:

Contents:

Main Unit, Main Unit Cover, Medication Container, Mesh Cap, Mask Adapter, Mouthpiece, 2 Batteries (AA/LR6) Adult Mask (PVC) Child Mask (PVC), Soft Pouch, Instruction Manual, Quick Setup Guide

#### **Value Added Features**

Delivery of Medication without Silicon Valve

Leaves Minimal Residual Medication

High Nebulization Rate

VVT Technology



# Value Added Features

2-in-1 compressor nebulizer and nasal aspirator

Alleviates congestion

Effective nebulization

Easy to use & clean





#### **Value Added Features**

Effective Nebulization Lightweight | User Friendly Easy to Clean | Low Noise







# **Wheeze SCAN**

WheezeScan uses **innovative listening technology** to identify wheeze quickly and accurately—so that parents have an in-the-moment tool to help them confirm what you've taught them to look and listen for.



**Uses advanced algorithm** to determine "wheeze" or "no wheeze" after being held below child's right collarbone for 30 seconds



**Clinically validated** to identify wheeze with 95% sensitivity and 90% specificity\*



**Helps parents remove uncertainty** anytime, anywhere in the moment of breathing distress and track events and triggers with free companion app

\*Clinically validated in trial of 600 children with suspected asthma or bronchitis in which panel of 2 to 3 doctors reviewed sound recordings to determine WheezeScan accuracy. Standard deviation of 0.219 sensitivity; 0.2948 specificity.




# **HOW TO USE NEBULIZER**





### ⚠ Caution:

Do not tilt the nebulizer kit at an angle of greater than 45 degrees in all directions. Medication may flow into the mouth or it may result in not nebulizing effectively.



# **IMPORTANCE OF PARTICLE SIZE**

ERS
European Respiratory Society recommends particle size (MMAD)
Recommendation as the most important parameter for selecting a nebulizer



## Why is MMAD important

Nebulizer creates an aerosol with a range of particle sizes and MMAD is the median diameter of the aerosol particle, which ensures efficient delivery of medication.

## **MESH NEBULIZER**

# 360° USABILITY

Unique 360° design allows for use in any position without spillage.





# **DIGITAL THERMOMETER**

#### Value Added Features

Fast 60 Sec Actual Measurement (Rectal)

**Accurate Measurement** 

Buzzer-Beeps when Measurement is Completed

Water-Resistant for Easy Cleaning





TECHNOLOGY BY OMRON HEALTHCARE IN JAPAN

#### Value Added Features

Quick Measurement-Oral, Rectal & Underarms

**Water-Resistant for Easy Cleaning** 

**Memory-Recalls Last Temperature** 

Buzzer-Beeps when Measurement is Completed

Flexible Tip



TECHNOLOGY BY OMRON HEALTHCARE IN JAPAN

#### Value Added Features

Quick Measurement in 1 sec

Gentle Non Intrusive Temperature Measurement

25 Memories Storage

LCDn Screen Display

2 mode-Silent & Buzzer-Beep

Selectable °C/°F mode

3-Type Measurement





INR 275/-

#### **Product Specifications**

Measurement Accuracy

± 0.1 °C [32.0 TO 42.0°C] ± 0.2°F [89.6 TO 107.6°F]

Measurement Range : 32.0 to 42.0°C (89.6 to 107.6 °F)

#### Contents

Test Battery Protective Case, Instruction Sheet, 8
Probe Covers





INR 399/-

#### **Product Specifications**

Measurement Accuracy

± 0.1 °C [32.0 TO 42.0°C] ± 0.2°F [89.6 TO 107.6°F]

Measurement Range: 32.0 to 42.0°C (89.6 to 107.6 °F)

#### Contents:

Test Battery (Alkline-Manganese Button Battery Lr41), Protective Case, Instruction Sheet, 5 Probe Covers

#### MC - 720-AP



**Medication Cup Capacity:** 

7 ml (cc) max

#### Measurement Accuracy:

Forehead measurement mode  $\pm 0.2^{\circ}$ C ( $\pm 0.4^{\circ}$ F) Surface measurement mode  $\pm 0.3^{\circ}$ C ( $\pm 0.5^{\circ}$ F) Room temperature mode  $\pm 2^{\circ}$ C ( $\pm 3.6^{\circ}$ F)

#### ontents:

Main unit, test battery (Cr2032 installed), probe lens



#### HBF - 214-AP



INR 4,860/-

#### **Product Specifications**

Weight accuracy:

2.0 kg to 40.0 kg: ± kg: ± 0.4 kg 40.0 kg to 150.0 kg : ± 1%

Power Source :

4 AAA Manganese batteries (R03) (You may also use AAA alkaline batteries (Lr03)

#### Package Contents:

Body composition monitor, 4 AAA manganese batteries (R03), instruction manual

#### **HBF - 224**





#### **Product Specifications**

Weight accuracy:

2.0 kg to 40.0 kg : ± kg: ± 400 g 40.1 kg to 135 kg : ± 1%

Power Source :

AAA Manganese batteries (R0)

#### Package Contents:

BCM, Trial battery C4 AAA (Manganes batteries) and Installation

#### **HBF - 222T-AP**



INR 5,950/-

#### **Product Specifications**

#### Weight accuracy

2.0 kg to 40.0 kg: ± kg: ± 0.4 kg 40.0 kg to 150.0 kg : ± 1%

#### Power Source :

4 AAA alkaline batteries (Lr03)

#### Package Contents:

Body Composition monitor, 4 AAA alkaline batteries (LR03), setup instructions, instruction manual, caution sheet

### **BODY COMPOSITION MONITOR**

### **Value Added Features**

Body Fat Percentage Visceral Fat Level Body Age Display Body Weight Skeletal Muscle % BMI (Body Mass Index)



Resting Metabolism

### Value Added Features

Body fat Percentage
Visceral Fat Level
Body Age Display
Body Weight
Skeletal Muscle %
BMI (Body Mass Index)
Resting Metabolism



#### **Value Added Features**

Body Fat Percentage
Visceral Fat Level
Body Weight
Skeletal Muscle %
BMI (Body Mass Index)
Resting Metabolism
Bluetooth Connectivity





#### Value Added Features

**Body Fat Percentage** 

Visceral Fat Level

Resting Metabolism

**Body Age & Weight Display** 

Segemental Skeletal Muscle Percentage

BMI (Body Mass Index)

Automatically device start measurement when you step on the device

**Auto Filtering Function** 

2 Sensors













#### HBF - 255T-APR



#### Product Specifications

2.0 to 100.0 kg with an increment of 0.05 kg 100 to 135.0 kg with an increment of 0.1kg

Power Source:

Package Contents:

HBF - 375-AP

Body composition monitor, 4 AAA alkaline batteries (LR03) Setup Instructions, Instruction Manual, Caution sheet Quick start Guide



**Body fat Percentage** 

Visceral Fat Level

Resting Metabolism

**Body Age & Weight Display** 

Segmental Subcutaneous Fat Percentage

Segmental Skeletal Muscle Percentage

4 User Profile

4 Sensors



#### TECHNOLOGY BY OMRON HEALTHCARE IN JAPAN



#### **Product Specifications**

0.0 kg to 40.0 kg : ±400 g 40.1 kg to 135.0 kg : ±1%

Power Source:

4 AA batteries (R6) (You may also use AA alkaline batteries (LR6)

**Package Contents:** 

Body composition monitor, 4 AA manganese batteries (R6). Instruction manual. Carry Bag

#### Value Added Features

**Body Fat Percentage** 

Visceral Fat Level

Resting Metabolism

**Body Age & Weight Display** 

Segmental Skeletal Muscle Percentage

BMI (Body Mass Index)

Automatically device start measurement when you step on the device

Segmental Subcutaneous Fat Percentage

- 4 Sensors
- 4 User Profile









### **HBF - 702T-AP**



#### Product Specifications

3.0 ti 100.0 kg with an incremer of 0.1 kg 100 to 135.0 kg with an incremer of 0.2 kg

Power Source :

4 A4 alkaline batteries (LR6)

**Package Contents:** 

Body composition monitor, 4 A4 alkaline batteries (LR6) Setup Instructions, Instruction Manual, Quick start Guide

# **WEIGHING SCALE**

#### **HN 289-EBK**



#### **Product Specifications**

Weight Range:

5 kg - 150 kg

Power Source :

CR2032 battery (Test battery is included)

Package Contents:

Digital Personal scale , One CR2032 battery for test use of the scale

# HN - 300T2-AP



#### **Product Specifications**

5 kg - 120 kg

Power Source:

4 AAA batteries

Package Contents:

Digital weight Scale, 4 "AAA" manganese batteries, Instruction Manual, Startup Guide

#### Value Added Features

Automatic On/Off

High Grade Tempered Glass Panel Sleek Ultra Thin Design **Precision Display** 













#### Value Added Features

Automatic ON/OFF Function 4 Sensor Accuracy Technology Precision Display to the closest 100a increment **Tempered Glass for Safety** 

**Auto Filtering Function** 

BMI on omron connect app



# **NERVE STIMULATOR**

#### Value Added Features

Compact

5 Region Massage

2 Massage Modes

Reusable Long Lasting pads

Frequency range Appx 1-108 Hz

Intensity level - 10



### Value Added Features

**Compact Design** 

Simple Operation

Reusable Long Lasting pads

3 Massage modes

6 Region Massage

Frequency range Appx 1-238 Hz

Intensity level - 15

**Carrying Pouch** 



#### Value Added Features

Frequency range Appx 1-1200 Hz

4 Massage MODES

Reusable Long Lasting pads

**Quick Relief button** 

Intensity level -10









OPTIONAL FOR ALL MODES Durable | Washable | Long Life



HV - F021

INR 4,990/-



OPTIONAL FOR ALL MODES Durable | Washable | Long Life



HV - F127

INR 5,990/-



OPTIONAL FOR ALL MODES Durable | Washable | Long Life



# **Electronic Nerve Stimulator Comparison Chart**

			•	
		HV-F013	HV-F021	HV-F127
*	Body Area Mode	5 Modes (Arms, Shoulders, Lower Back, Legs, Feet)	6 Modes (Arms, Shoulders, Lower Back, Joints, Legs, Feet)	4 Modes (Shoulders, Lower Back, Legs, Feet)
	Massage Mode	STEADY KNEAD	TAP RUB	TAP RUB  KNEAD PUSH
<b>9</b>	Frequency Range	Approx. 1 to 108 Hz	Approx. 1 to 238 Hz	Approx. 1 to 1200 Hz

# **OXYGEN CONCENTRATOR**

# Oxygen Concentrator

Efficient and Stable

90-96%

High Oxygen Concentratration



#### 2 in 1 functionality

Oxygen Production







INR 79,990/-

#### Value Added Features



Stable oxygen supple 24/7

≥90%

High Oxygen Concentration

### 4 Major Safety Alarms

- - Faulty alarm
- Low oxygen alarm
  - ) Power failure alarm
  - $\overset{\smile}{\mathbb{T}}$  o
    - Overheating alarm



# **Timed Oxygen Production Mode**



Set the oxygen inhalation time as you like



You can set it flexibly according to your own needs



Up to 8 hours can be set



You can view the single/ cumulative use time at the same time



**Advanced Features** 

## **PROFESSIONAL CARE**

#### Value Added Features

5 Cuff Sizes Available

**Dual Measurement Mode** (Oscillation + Auscultation )

Easy Maintenance

Irregular Pulse indicator

**Motion Stop Function** 



CLINICALLY VALIDATED



#### Value Added Features

Dual Measurement Mode (Oscillation + Auscultation)

**Ergonomic Design** 

5 Cuff pieces

Power Saving auto power off

Backlight LCD

Last Reading Display

Rechargeable battery



CLINICALLY VALIDATED



#### Value Added Features

Moveable arm Cuff

Accurate measurement navigation

Elbow detection sensor

Cuff Range 17-42cm

Various print out formats

**LED Display** 

**USB Data transfer** 



CLINICALLY VALIDATED



#### HBP - 1120-IN



# **Product Specifications**

Pressure 0 to 299 mmHg Pulse 40 to 180 beasts/min

Memory

Pressure ±3mmHg Pulse ±5% 0f display reading

#### HBP - 1320-AP



#### Product Specifications

Measurement Range:

Pressure 0 to 299 mmHg Pulse 40 to 180 beasts/min

Memory

**HBP - 9030** 



INR 1,69,490/-

#### Product Specifications

Measurement Range:

Pressure 0 to 299 mmHg Pulse 40 to 180 beasts/min

Memory

# **PROFESSIONAL CARE**

#### **NE-C900-AP**



#### **Product Specifications**

weight accuracy :

2.0 to 100.0 kg with an increment of 0.05 kg 100 to 135.0 kg with an increment of 0.1kg

Particle Size

3 µII

#### Package Contents:

Body composition monitor, 4 AAA alkaline batteries (LR03) Setup Instructions, Instruction Manual, Caution sheet Quick start Guide

# Value Added Features

Ergonomic design
Effective Nebulization delivery
Appropriate Particle size
Compressor based power
performance



Easy to clean

### **NE - U780-AP**



#### **Product Specifications**

**Medication Cup Capacity:** 

10ml minimum and 150ml maximum

Particle Size

4.5 µm

#### Package Contents:

Main Unit, Removable water tank, Medication cup cover, Medication cup holder, Inhalation hose M, Mouthpiece, Medication Cup [2 Pcs] Power Cord, Instruction Manual.

#### Value Added Features

Efficient medication delivery

Nebulization rate and

Adjustable Airflow volume

Effective Nebulization

Appropriate Particle size

Powerful Machine

Removable Water tank

Audible Alarm

madible marin

Large Display

**Easy Cleaning** 



#### Accessories For Blood Pressure Monitor

#### **AC ADAPTER S**



Compatible with all Blood Pressure Monitor models In case you don't want to use batteries, purchase your AC Adaptor for continuous use.

Model Code - HHP-BH01-IN

Medium Cuff (fits arm circumference 22-32cm) Model Code - HEM-CR24-BAP



#### **SMALL CUFF Cs24**



Small Arm cuff for Home Bp Monitor that's compatible with devices except HEM 7600T, HEM 7361 T, and HEM 7156/A.

This arm cuff can comfortably fit on the arm sizes ranging between 17 to 22 cms. Model Code - HEM-CS24-BAP

Compatible with only Blood Pressure Monitor Models HEM 7361-T & HEM 7156/A Model Code - HEM-FL31-BAP



### WIDE RANGE CUFF (M-L)



INR 890/-

Compatible with all Home BP Monitor Model's Except HEM 7600T/ HEM 7361-T & HEM 7156/A Arm Cuff for BP Monitors.

Arm Size 22 cm to 42 cm Model Code -HEM-RML31-BAP

Accessories For Blood Pressure

### **NEBULIZER AIR TUBE** (NEC-28)



Compatible Model No.: NEC-28 Nebulizer Air Tube

Model Code :9956270-0

Compatible Model No.: NE-C101 Nebulizer Kit, Child Mask & Air Tube Set

Model code: NEB-NSET5-11AP



#### NEB-NSET5 - (11AP)



Compatible Model No.: NE-C101 Nebulizer Kit. Adult mask & Air Tube Set

Model code: NEB-NSET4-11AP

Compatible Model No.: NE-C28/C29, NE-C30 Set of Mouthpiece & Med Cap

Model Code - 9956268-9

# WIDE RANGE CUFF (M-L)

INR 250/-

### NEBULIZER KIT SET -(NEC-803)



Compatible Model No.: NEC-803 Nebulizer Kit

Model code: 9956268-9

Compatible Model No.: NEC-106 Nebulizer Kit

Contents: Nebulizer Kit, Air tube, Mouthpiece, Adult mask, Child Mask, Air filterx5

Model code: NEB-ASKIT-16AP

**NEBULIZER KIT SET** (NEC-106)



INR 449/-

Accessories For Nebulizers

#### **ADULT MASK PVC** (NEB-MSLP-E)



Compatible Model No.: NE-C28/C29, NE-

C30 Adult Mask Only

Model Code: 9956275-1

Compatible Model No.: NE-C28/C29, NE-

C30 Child Mask Only

Model Code :9956276-0

**CHILD MASK PVC-**(NEB-MSMP-E)

INR 150/-

#### **AIR FILTER (X5) - (C28-2-C1)**



Compatible Model No.: NE-C28/C29, NE-C30 Air Filters for nebulizers

Model Code: 9956271-9

Compatible Model No.: NE-C101 Air

Filters for nebulizers

Model Code: NEB-AFR-30AP

AIR FILTER (X5) - (NEC101)



### AIR FILTER (X5) - (C106)



Compatible Model No.: NE-C106 Air

Filters for nebulizers

Model Code: NEB-AFR-30AP

Replacement Aspirator Accessory Set for the OMRON DuoBaby Nebuliser.

Model code: NEB7000

NE-C301 (Duo baby) **Aspirator Accessory** 



### Accessories For Oxygen Concentrator

### Y-506W (HAO-5000M-CAAP)



Nasal Oxygen Canula

Model code: MB01-L-2.0

#### Y-506W (HAO-5000M-NKAP)



INR 545/-

Nebuliser mask for single use

Model code: MW01-L-2.0

## **Accessories For Tens Pain Therapy Devices**

Compatible Model No.: HV-F013, HV-F021, HV-F127, HV-F128Pad, 1 Pair for replenishment in case the stickiness of pads is reduced after continuous use with Nerve Stimulator HV-F013

Model Code: HV-LLPAD-C1

LONG LIFE PAD (1 PAIR)



Accessories For Professional Care

# NEBULIZER KIT SET WITH ADULT MASK (NEB-NSTE4-90AP)



NE-C900 Nebulizer Kit, Adult Mask & Air Tube Set

Model code: 9063553-5

# NEBULIZER KIT SET WITH ADULT MASK (NEB-NSTE4-90AP)



NE-C900 Nebulizer Kit, Child Mask & Air Tube Set

Model code: 9063554-3

# NE-C301 (Duo baby) Nebulizer Kit



Replacement Nebulizer Accessory Set for the OMRON DuoBaby Nebuliser.

Model Code: NEB2012

## OMRON CUSTOMER SUPPORT LOCATION (WEST) Experience and Service Centre

#### Thane

Bafna Surgicals, 3, Landmark Arcade, Besides Axis Bank, Louiswadi Service Road, Thane (W), Pin Code - 400604 M: +91 9555795795

#### Ahmedabad

Infinity Healthcare A-101, SUMEL-6, Jupiter Mill Compound, Opp. Hanumanpura, Shahibaug Dudheshwar Road, Ahmedabad, Gujarat, 380004 t-9018642027, 9840821663

#### Mumbai

ACCEL IT SERVICES Unit no. 821, 2nd Floor, Solitare Corporate Park Andheri - Ghatkopar Link Road, Chakala, Andheri (East), Mumbai - 400093 T - 022 62040901

#### Vadodara

Arihant surgical co FF/106 benison complex opp rajlaxmi society near chakli circle oproadvadodara-390007 mob-7016106315

#### Pickup Centres

#### Ahmedabad

ACCEL IT SERVICES 201/B, City Pride Near Hotel Nalanda, Mithakhali Cross Road Navarangpura Ahmedabad, Gujarat Pin Code - 380006 +91 79 26561927/28. +91 7824819590

#### Pune

ACCEL IT SERVICES
Office No 2, Ground floor,
Mayfair Towers 2, Wakadewadi,
Shivajinagar, Pune,
MAHARASHTRA-411005,
T- 020-41204835

#### Jabalpur

Varsha Surgical & Medicals BK Tower, Old Savitri Talkies Compound Ukhri Road, Baldleo Bagh, Jabalpur Madhya Pradesh Pin Code - 482002 +91 900 9437437

#### Indore

ACCEL IT SERVICES Ground Floor, Malwa Tower 10, Block # 5, Old Palasia, Guitarwala Square, A.B.Road, Indore, Madhya Pradesh Pin Code - 452001 +91 731-4995347

#### Nagpur

ACCEL IT SERVICES Plot No 49, Bhagva Ghar, Beside Traffic Children's Park, Near Tadka Restaurant, Dharampeth, Nagpur Pin Code - 440010 +91 712-2995131

#### Pune

Ronak Pharma Shop no-53&54, Mayur Center, Opp. Unique Hospital, Dawa Bazar, Chinchwad, Pune-411019 +91 9075894107

#### Goa

Drogaria Raicer Sunshine Bldg, Jose Inacio Loyola Rd, Govind Poy Raiturcar, Margao,GOA - 403601 T-9890630763

#### RAIPUR

Shree Laxmi Surgical Near Lotus Showroom , Opp. Hotel Arraz, Raipur-492001 T-8889486024

#### Pune

Nitin Sales Shop No-1+ Godown, Shree Dattakrupa Society Soc. H.No.1411,Sadashiv Peth, Pune-411030 T - 7774932291

#### Mumbai

Khimasia Drugs & Surgical House G/86-A Gosar Niwas Behind Shah Jayantilal Gosar & Co. Juhu Church Road Juhu, Mumbai Pin code: 400049 +91 9930591294 Surat ACCEL IT SERVICES

ACCEL IT SERVICES
Office No.M/20, 21, 22, 23
Shreenath Complex
Opp. Jivan Jyot Honda Showroom
Udhna Sachin Road, Surat, Gujarat
Pin Code - 394210 +91 261 2272571/72

#### Baroda

ACCEL IT SERVICES Ltd 001, Yogiraj Apartment, 45 B Nutan Bharat Soc, Nr. M. K. High School Alkapuri, Baroda, Gujarat Pin Code - 390007 +91 9033050153/ 9033060155

#### Nasik

ACCEL IT SERVICES Plot No: 33, Ground Floor, Behind RBS Bank, Krushi Nagar College Road, Nasik, Maharashtra Pin Code - 422005 +91 253 3023103, +91 253 2310177

#### Bhopal

Smith System Plot-194 (G216-217), Surya Towers, (Beside Sahara Kunj & Opposite DB Mall). MP Nagar Zone-1 Bhopal (MP), PIN CODE - 462011 T - 0755-429554 ,9893720720

#### Nagpur

Chandra Enterprises 104, Guru Balwant Complex,1st Floor, Near Yeole Tea Stall, Lokmat Sq., Nagpur-440012 T-9552146982/0712-2427748

# OMRON CUSTOMER SUPPORT LOCATION (EAST) Service Centre

#### Kolkata

ACCEL IT SERVICE P-89, CIT Road Entally, Puddapukur, Kolkata- 700014 Ph: +91 33 40325607, +9133 40325615

#### **Pickup Centres**

#### Bhubaneshwar

ACCEL IT SERVICES 1st Floor, Plot no. 100, Saheed Nagar, Dist: Khurdha Bhubaneshwar, Odisha Pin Code - 751007 91 674 2548540. +91 674 2548492

#### Howrah

Air Voice 374, Belilious Road, Opp: Sales Tax Building Near Maidan, Howrah, West Bengal Pin Code - 781007 +91 9748347344, +91 9230215888

#### Patna

ACCEL IT SERVICES 2nd Floor, Flat No. 2/5, Surya Apartment, Fraser Road, PATNA, BIHAR-800001 , T-0612 -2216083

#### Ranchi

Akshat agencies Indrapuri Road No-03,Ratu Road, Ranchi-834005 Mob no-8340266616, 9386318567, 9431169638.

#### Rourkela

Om Sai Meditech Plot no-273/530, Near SBI Bank, 21-Shakti Nagar, Rourkela -769014 Dist -Sundargarh, Odisha M: +91 7327820309 / 9583422620

#### Siliguri

Surgical Hub Pradhan Nagar Meghnad Sarani, Siliguri T: 0353 35364199

#### Jamshedpur

Bajwa Surgicals H.No. 2, Line No. 1, A-Block, Tuiladungri ,Golmuri , Near R.D. Tata Golchakkar, Jamshedpur, Jharkhand - 831003 +91 7209046879

#### Siliguri

Sarovar Surgi Co. Sisir Bhaduri, Khudirampally, Opp. Kanchenjunga Stadium, Near Badshah Biryani, Siliguri, West Bengal Pin Code - 734001 +91 9830780009

#### Guwahati

ACCEL IT SERVICES Mitra Building, Ashram Road Ulubari Chariali, Guwahati, Assam Pin Code - 781007 +91 361 2450611/12

#### Raipur

Shree Laxmi Surgicals Shop No-6, Main Road medical complex rajbandhan maidan raipur-492001 M: +91 8889486024

#### Bilaspur

SHRI BALAJI SERVICES C/o Shri Balaji Surgical, Medical Complex, Telipara, Bilaspur, Chattisgarh-495001 M: +91 9755559018

# OMRON CUSTOMER SUPPORT LOCATION (SOUTH) Experience and Service Centre

#### Chennai

Accel IT Services (A division of Accel Limited), 1st Floor, SFI Building, No. 175, Valluvarkottam High Road, Chennai - 600034, T +91 7824845801

#### Chennai

The Health And Wellness Shoppe No 55 Hunters road Choolai - 600112 Landmark - Near Doveton Bus Stop, Opp To St Paul Church Landline - 044-79662925, Mobile - 9080226851

#### Hyderabad

New Prince Agencies , Plot no-5/1 & 88 , Chaitanya Chambers , Chaitanyapuri ,Hyderabad-500060. T-040 4953081, M-9949945335

#### Cochin

ACCEL IT SERVICES Door No-38/261 Behind Eastern Corporate Office, Edappally Arakkakadavu Road, Edappally P.O. Edappally, Kochi, Kerala 682024. T- 0484 4863201/0484-4863203

#### Bangalore

Jaij Marketing 109/6, 2nd floor,20th main road, West of Chord road, 5th block, rajaji naga, Bangalore-560010 Next to Variar bakery Mobile No – 09148331117

### **Pickup Centres**

#### Hubli

ACCEL IT SERVICES Ground Floor, Diamond Corner Opp Sawai Gandhrava Hall Deshpande Nagar, Hubli - Karnataka Pin Code - 580029 +91 836 2213755, +91 836 4256517

#### Kannur

Ecom Surgicals CW35/1672, 3rd Floor Jubilee Bazar, Onden Road, Kannur-670001 M: +91 9447046636

#### Mysore

Mr Majula Prakash #1079, Railway station road,chamrajpura, Mysore-570005 M: +91 9900478385

#### Gulbarga

Mr Syed , Bl complex, opp. Gandhi nagar, Kaman Humnabad road, Gulbarga-585101, T-9845414148

#### Mangalore

HCiT Services 15-18, "Astral", Padavu, Kalpane Kulashekar, Mangalore- 575005 M: +91 9686676755

#### Kollam

Spartan impex KP10/124E BYPASS ROAD Opp Malayalam motors (skoda showroom) Palkulangara.kollam-691004 M: +91 9544000786

#### Coimbatore

S.k.healthcare, Door No.113-a, Front Portion, Balaji Nagar, New Siddhapudur, Coimbatore, Tamil Nadu-641044, T- 0422-4373668, 9380669668

#### Shivmoga

Mr Umaravi shankar #9, Salur matta apartment, 1st cross basavanagudi, opposite Vinayaka park, shimoga-577201, M: +91 9686516276

#### Bangalore

Jaij Marketting & Services No.-947/948, ITPL Road ,rajaplaya, Opp. Punjab National Bank, Hoody Village, Banglore-560048 , T - +91 80-41200604, +91 9880411117

#### Coimbatore

ACCEL IT SERVICES No-94A & 95, Prasad Chambers East Arokiaswamy Road, RS Puram Coimbatore - Tamil Nadu Pin Code - 641002 +91 422 4976306

#### Visakhapatanam

ACCEL IT SERVICES Rednam Regency 47-10-10, 3rd Floor,T2 Dwarakanagar 2nd lane, Visakhapatnam-530016 T + 91 891-2567881

#### Puducherry

ACCEL IT SERVICES No 12, 1st Floor, 2nd Cross Road Jawahar Nagar, Boomianpet, Puducherry Pin Code - 605005 +91 413 2200314, +91 413 2200315

#### Salem

ACCEL IT SERVICES 37/6C, Ramakrishna Road Shanmugam Colony, Salem, Tamil Nadu, Pin Code - 636007 +91 427 4300073, +91 427 2336457

#### Panaji-Goa

ACCEL IT SERVICE S No: 501/521, 5th Floor, Citi Centre Patto Plaza, Panjim, Panaji Goa, Pin Code - 403001 +91 832 6634501, +91 832 6634502

#### Trichy

ACCEL IT SERVICE S C-87, North East Extension, 6th Cross, Sastri Road, 1 st Floor, Nagappa Complex Thillai Nagar, Trichy Tamil Nadu, Pin-620018 +91 431 4051630, +91 431 4011333

#### Thiruvananthapuram

ACCEL IT SERVICE S TC-15/1805(3), I Floor, Mahesh Estates Vazhuthacaud, Thiruvananthapuram Kerala, Pin Code - 695014 +91 471 4014494

# OMRON CUSTOMER SUPPORT LOCATION (NORTH) Experience and Service Centre

#### Gurugram

MS Services Opp Krishna Mangalam Marrige garden Sheetla Mata Road, Gurugram Pin Code -122001

Tel: +91 7048999792, +91 9212733632

#### Ludhiana

Vaishnavi Surgical Shop no.-26, 1st Floor, JMD mall near Railway station, Ludhiana-141008, M: +91 9872857957

#### Varanasi

Balwant Pharma C-27/217, JAGATGANJ, Varanasi -221001, M: +91 9935507271, 7897271786

#### Mohali

Office Plus limited SCO-12 Sector-68 Mohali, Near Punjab & Sind bank Punjab-160062 T-0172 4784666, 9888200675

#### New Delhi

ACCEL IT SERVICES A-61/1, Ground Floor, Okhla Industrial Area, Phase-II, New Delhi, Pin Code - 110020 Tel: +91 11 46575100, +91 11 46575112

#### Lucknow

Docket Care systems 52,53 Chakarpuri Colony Paper Mill Road, Nishatganj, Lucknow-226006 M: +91 9415008355, 9918135014

#### Jaipur

Raj Associates B-61-B,Ved Vila Colony-B,Swej Farm, Ram Nagar,Sodala Jaipur-302019 M: +91 7014884010, 9982847772

### **Pickup Centres**

#### Agra

Vaishno Pharmaceutical 37/264,Opp. Vijay market, Dayal Bagh Road New agra ,Agra 28200, T -9917018877

#### Ghaziabad

VYAS ENTERPRISES , Shop No-6, Sec-10, D/155, Kalpana Complex,Vasundhara, Ghaziabad-201012 , T: +91 120-4142206, M:+91 9873169612

#### New Delhi

Omkar Services first floor, U-197, Shakarpur, behind Laxmi nagar Bus Stand in front of metro pillar no-32, delhi-110092 M: +91 7053168983, 9311831923

#### New Delhi

Omkar Services Plot No.140, First Floor, Khasra No- 9/5 Kakrola Housing Complex, Dwarka Delhi, Pin Code - 110078 +91 7838326175, +91 7291938617

#### New Delhi

Medicare Surgical 330, Patparganj Indl. Area, Near Ola Service Center, Delhi-110092 T 8077413922

#### Faridabad

National Services 5B/8BP, Railway Road NIT Faridabad Near Banke Bihari Mandir, 121001 +91 1294311261, 9319965879

#### Amritsar

Konica Minolta Authorised service center (Printer) Malviya Complex, Malviya Road, INA Colony ,Near Sewa samiti hospital, Amritsar, Punjab-143001 M: +91 9815503055

#### Dehradun

ACCEL IT SERVICES No-109 Tagore Villa Mauza Chukkuwala Mehar Colony, Chakrata Road Dehradun- Uttrakhand Pin Code - 248 001 +91 135 2655797

#### Hisar

S.L Services 519B, Krishna Nagar, Camp Chowk, Near Prem Bandhan Banquet Hall, Hisar - 125001 +91 9254289838, +91 9034770099

#### Ambala

Tandon Surgical Shop No. 5438/3,cross Road, No. 4, Nicholson Road, Ambala Cantt. - 133001 T-9812144534/9813204726/01712654424

#### New Delhi

OMKAR SERVICES 487/7, Peera Garhi Chowk, New Delhi -110087 T- 9310126767

#### Kota

Sharma Electronics 384 Shopping Centre, Near IDBI Bank Kota, Rajasthan Pin Code - 324006 M: + 91 8947037114/9829116833

#### Jodhpur

Vijay Communication 103, Behind New Kohinoor Cinema Masuriya, Jodhpur, Rajasthan Pin Code - 342001 +91 7568705605, +91 9928715825

#### Kanpur

ACCEL IT SERVICES 112/8, Benajhabar Main Road, Opposite Sachiv Nivas, Kanpur - Uttar Pradesh Pin Code - 208002 +91 8939822084

#### Chandigarh

Accel IT Services SCO Number-5, 1st floor, Sector-20D, Chandigarh-160020 Tel No-0172-3023003

#### Srinagar J&K

Swift Orthotics Near JLNM Hospital, Raina wari , Srinagar -190003 T-0194 2507152 , 7006189139

#### Sirsa

BANSAL PHARMACEUTICALS, OPP. R. S.D. SCHOOL, MONGA DENTAL WALI GALI, DABWALI ROAD, SIRSA-125055 T-9812600629/9992022529/8901119329

#### Kurukshetra

SARASWATI BIO CHEM , opp. Gita Sen . sec. school, RLY Road, Kurukshetra -136118 T-9812700403

#### **GHAZIABAD**

Radiant Corp Shop No-lg-3, Lower Ground Floor, Plot No-c 17, Pandit Plaza, Sector-13, Vasundhara, Ghaziabad-201012 T-9561296397

#### Noida

Medicare Solution G 3 Gopal Ji Barat Ghar Village Shahpur, Near Jaypee Hospital, Sector 128, Noida, Uttar Pradesh Pin Code - 201304 +91 931 0023407

#### Delhi

Medicare Surgicals 520 Bhagirath Palace, Near Central Bank of India, Jaskan Pal Building, Delhi Pin Code - 110006 +91 8077413922

#### Udaipur

Unique Medicos Shop No. 4, Ground Floor, Shrinath Plaza, Opp. M.B. Govt. Hospital, Chetak Cirlcle, Udaipur - 313001 02942417416, +91 9352526141

#### Rohtak

AROMA SURGICAL HOUSE, Delhi road, model town, Rohtak-124001. T-9812044476/9896066619

#### Sonepat

Saatvik Medicos 32-A Geeta Bhawan Sonepat-131001 T-9812190008 / 9728100000

#### Rewari

YADAV SURGICAL 299/23, Bharawas Road, Near Dr Atamprkash Hospital Company Bagh, Rewari -123401 T-9413463632

#### New Delhi

Decent Healthcare 53/4, Ashok Nagar, Near Swastik Hospital, New Delhi-110018 T-9136230343

#### Jammu

Jogindra pharmaceuticals 151 raghunathpura, jammu 180001 +911912573357

# For more Product Information or Demo, Kindly Contact undersignedn:

Name
Phone No
E-mail Id
Location/City



### **OMRON HEALTHCARE INDIA PVT. LTD.**

6th Floor, B-Block, Sewa Tower, Plot No. 19, Sector 18, Udyog Vihar, Maruti Industrial Complex, Gurugram, Haryana-122008

Toll Free: 1800-419-0492

Monday to Friday-9 am to 5:30 pm (Except Holidays)

Email: callcentre@omron.com
Website: www.omronhealthcare-ap/in | www.omronbrandshop.com



